




# Duckworth & Kent Ltd.


titanium surgical instrument manufacturer

## INSTRUMENT CARE & HANDLING



Your Instrument Care  
is Patient Care





ALWAYS REMEMBER TO CHECK FOR



Irregular/Bent tip




Chipped blade





Tip bent & damaged




Missing graduation



Improperly aligned tip





Irregular/Bent tip




Irregular/Bent tip





Irregular/Missing point



Irregular/Chipped tip



Irregular/Bent tip



Irregular/Bent tip



Irregular/Bent marker