



Duckworth & Kent

Titanium Reusable Ophthalmic Instrument Manufacturer



A CLASSIC DESIGN

REDEFINED




The **Heritage Collection** is a set of four highly versatile forceps that have been re-engineered with the latest innovations in single piece titanium construction, first introduced by our E Range Forceps.

Offering the same **durability, stability** and **exceptional quality** that our customers have come to expect from our precision-engineered forceps, this concise collection provides a greater choice by reintroducing a classic popular handle design.



Giving You More Choice Handle Styles

E Range Handle




Heritage Handle




2-100 **Pierce Notched Forceps**

£155.00




- Pierce 0.25mm notched, 6mm tying platforms
- Straight shafts
- Short flat handle, overall length 89mm
- Single piece instrument





2-110 **Bonn Toothed Forceps**

£155.00




- 0.12mm, 1 into 2 teeth, 6mm tying platforms
- Straight shafts
- Short flat handle, overall length 89mm
- Single piece instrument




2-132 **Troutman-Barraquer Toothed Forceps**

£185.00




- 0.12mm, 1 into 2 teeth, 6mm tying platforms
- 2mm colibri shafts
- Short flat handle, overall length 89mm
- Single piece instrument



2-500 **Tying Forceps**

£149.00



- 6mm tying platforms
- Straight shafts
- Short flat handle, overall length 89mm
- Single piece instrument

